



HYPHEN
BioMed

155 rue d'Eragny, 95000 Neuville-sur-Oise, France

Tél : +33 (0)1 34 40 65 10

Fax : +33 (0)1 34 48 72 36

www.hyphen-biomed.com

info@hyphen-biomed.com



BIOPHEN™ Plasma Calibrator

REF 222101

CAL 12 vials x 1 mL

IVD

English, revision: 11-2025

INTENDED USE:

For calibration of hemostasis assays using a quantitative automated method.
This device of *in vitro* diagnostic use is intended for professional use in the laboratory.

SUMMARY AND EXPLANATION:

Technical:

The following table shows the various parameters, which are measured using assays from HYPHEN BioMed, and according to the package inserts:

Parameters	Method
Fibrinogen (Fbg)	Activity / Antigen
$\alpha 2$ – Antiplasmin ($\alpha 2$ -AP)	Activity
Antithrombin (AT)	Activity / Antigen
Factor II (FII)	Activity
Factor V (FV)	Activity
Factor VII (FVII)	Activity
Factor VIII (FVIII)	Activity
Factor IX (FIX)	Activity
Factor X (FX)	Activity
Factor XI (FXI)	Activity
Factor XII (FXII)	Activity
Factor XIII (FXIII)	Activity
Plasminogen (Plg)	Activity
Protein C (PC)	Clotting activity
Protein C (PC)	Chromogenic activity
Protein C (PC)	Antigen
Protein S (PS)	Clotting activity
Free Protein S (Free PS)	Antigen
von Willebrand Factor Antigen (vWF: Ag)	Antigen

REAGENTS:

CAL Normal citrated human plasma, lyophilized.

Calibrators contain stabilizing agents.

The calibrator concentrations may vary slightly from one batch to another. For the assay, see the exact values indicated on the flyer provided with the kit used.

The product is classified as non-hazardous and is not subject to labeling according to EC Regulation No. 1272/2008 [CLP].

WARNINGS AND PRECAUTIONS:

- Some reagents provided in these kits contain materials of human origin. Whenever human plasma is required for the preparation of these reagents, approved methods are used to test the plasma for the antibodies to HIV 1, HIV 2 and HCV, and for hepatitis B surface antigen, and results are found to be negative. However, no test method can offer complete assurance that infectious agents are absent. Therefore, users of reagents of these types must exercise extreme care in full compliance with safety precautions in the manipulation of these biological materials as if they were infectious.
- Waste should be disposed of in accordance with applicable local regulations.
- Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.
- Summary of Safety and Performance (SSP) is available in the European database on medical devices (see Eudamed public website: <https://ec.europa.eu/tools/eudamed> or on request to HYPHEN BioMed).

REAGENT PREPARATION:

Gently remove the freeze-drying stopper, to avoid any product loss when opening the vial.

CAL Reconstitute the contents of each vial with exactly 1 mL of distilled water.

Shake vigorously until complete dissolution while avoiding formation of foam and load it directly on the analyzer following Application Guide instruction.

This plasmatic reagent can be more or less turbid after reconstitution. This turbidity is mainly due to plasma lipids that, after freeze-drying, become "less" soluble and may form a slight deposit. If necessary, let each vial stabilize 10 minutes at room temperature and shake before use.

STORAGE AND STABILITY:

Unopened reagents should be stored at 2-8°C in their original packaging. Under these conditions, they can be used until the expiry date printed on the kit.

CAL Reagent stability after reconstitution, free from any contamination or evaporation, and stored closed, is of:

- For AT, PC, FVII+X, Plasminogen, FII, FVII, FIX, FX, FXI, FXII, FXIII, Fibrinogen, $\alpha 2$ -AP, Free PS and vWF: Ag:
 - 24 hours at 2-8°C.
 - 28 days frozen at -20°C or less*
 - Stability on board of the analyzer: see the specific Application Guide.
- For FVIII, FV and PS:
 - 8 hours at 2-8°C.
 - 28 days frozen at -20°C or less*

- Stability on board of the analyzer: see the specific Application Guide.

*Thaw only once, as rapidly as possible at 37°C and use immediately.

REAGENTS AND MATERIALS REQUIRED BUT NOT PROVIDED:

- Laboratory material.
- HYPHEN BioMed hemostasis reagents.
- CS-series, CN-series, CA-series*, STA-R® family, ACL-TOP® family.
- * only when used with reagents claiming CA series.

TRACEABILITY:

Lot to lot variability measured on 3 lots is: %CV $\leq$ 10%.

Calibrator is traceable to WHO International Standards in force for Fibrinogen, Antithrombin, FII, FV, FVII, FVIII, FIX, FX, FXI, FXII, FXIII, Protein C, Protein S, Free Protein S and von Willebrand Factor parameters.

Calibrator is traceable to normal plasma pool for $\alpha 2$ -AP and Plasminogen parameters.
Certificate of traceability and uncertainty is available on the HYPHEN BioMed website:

Parameters	Uncertainty	Parameters	Uncertainty
Fibrinogen (Fbg)	± 0.23 g/L (antigen) ± 0.17 g/L (activity)	Factor X (FX)	$\pm 6\%$
$\alpha 2$ – Antiplasmin ($\alpha 2$ -AP)	$\pm 9\%$	Factor XI (FXI)	$\pm 10\%$
Antithrombin (AT)	$\pm 5\%$ (antigen) $\pm 8\%$ (activity)	Factor XII (FXII)	$\pm 11\%$
Factor II (FII)	$\pm 4\%$	Factor XIII (FXIII)	$\pm 9\%$
Factor V (FV)	$\pm 9\%$	Plasminogen (Plg)	$\pm 5\%$
Factor VII (FVII)	$\pm 6\%$	Protein C (PC)	$\pm 4\%$ (antigen) $\pm 9\%$ (activity)
Factor VIII (FVIII)	$\pm 9\%$	Protein S (PS)	$\pm 10\%$
Factor IX (FIX)	$\pm 9\%$	Free Protein S (Free PS)	$\pm 5\%$
		von Willebrand Factor Antigen (vWF: Ag)	$\pm 6\%$

QUALITY CONTROL:

For calibration of hemostasis assays.

The target values are determined from multi-reagent and multi-instrument tests.

The use of quality controls serves to validate method compliance, along with between-series assay homogeneity for a given batch of reagents.

Include the quality controls with each series, as per good laboratory practice, in order to validate the test.

A new calibration curve should be established, preferably for each test series, and at least for each new reagent batch, or after analyzer maintenance, or when the measured quality control values fall outside the acceptance range for the method.

LIMITATIONS:

- If the calibrators are used under measurement conditions other than those validated by HYPHEN BioMed, the test results may vary. The laboratory is responsible for validating the use of these calibrators in its own analytical system.
- Any reagent presenting no limpid appearance or showing signs of contamination must be rejected.

e-IFU (other languages) are available on www.hyphen-biomed.com.

For customer support and Application Guides, please contact your local provider or distributor (see www.hyphen-biomed.com).

Changes compared to the previous version.

The following symbols may appear on the product labeling:

REF	Catalogue number	LOT	Batch code	IVD	In-vitro diagnostic medical device
Rx	Numerical < x> identification of reagent	i	See instructions for use	WHO STD	WHO standard code
Temperature limitation		Manufacturer		YYYY-MM-DD	Use by
CE XXXX	CE marking of conformity with notified body ID number.	→	Reconstitution volume	CONTENTS	Contents
Cx	Numerical < x> identification of control	i-MA	See instructions in Method Application guide	CONTAINS	Contains
EXP	Expiration date	Σ	Contains sufficient for <n> tests	UNIT	Measurement unit
TARGET VALUE	Target Value	Keep away from sunlight and heat		CALx	Numerical < x> identification of calibrator
UDI	Unique Device Identifier	BIO	Contains biological material of animal origin		Contains human blood or plasma derivatives
DANGER	Danger	WARNING	Warning	UK CA	UKCA marking of conformity
Biological risks		ACCEPTANCE RANGE	Acceptance range		



HYPHEN BioMed

155 rue d'Eragny, 95000 Neuville-sur-Oise, France

D750-02/BI/2101/v6