Sekisui Diagnostics, LLC
500 West Avenue, Stamford, CT 06902
Tel. (203) 602-7777 Fax (203) 602-2221

LAtrol®

REF 816A/816N

american diagnostica GmbH
Kaplaneigasse 35, Pfungstadt, D-64319 Germany

REVISION B

REF 816A: 10 vials each containing 0.5 mL of LATrol™ Abnormal Control
10 Fläschchen mit je 0.5 mL LATrol™ Abnormale Kontrolle
10 flacons contenant chacun 0,5 mL de LATrol™ Contrôle Anormal
10 viales con 0,5 mL de Control Anormal LATrol™ cada uno
10 flaconi contenenti ciascuno 0,5 mL di LATrol™ Controllo Anormale

REF 816N: 10 vials each containing 1.0 mL LATrol™ Normal Control
10 Fläschchen mit je 1.0 mL LATrol™ Normal Kontrolle
10 flacons contenant chacun 1,0 mL de LATrol™ Contrôle Normal
10 viales con 1,0 mL de Control Normal LATrol™ cada uno
10 flaconi contenenti ciascuno 1,0 mL di LATrol™ Controllo Normale

ENGLISH

INTENDED USE

The LATrol™ Abnormal Control (REF 816A) and LATrol™ Normal Control (REF 816N) plasmas have been developed for use as part of daily quality control procedures for Lupus Anticoagulant (LA) testing. These control plasmas are designed to be used with ACTICLOT® dPT™ (REF 824), DVVtest® (REF 810/825), DVVconfirm® (REF 815/815L) and ACTICLOT® LA Complete™ (REF 826).

DESCRIPTION

LATrol Abnormal Control is a lyophilised preparation of a Lupus Anticoagulant plasma with added buffer, which has been determined to be positive for LA in accordance with the revised criteria of the Scientific and Standardization Committee's (SSC) Subcommittee for the Standardization of Lupus Anticoagulants.¹

LATrol Normal Control is a lyophilised preparation of a multi-donor normal plasma pool with added buffer. This normal plasma may also be used in mixing tests to determine the presence of inhibitors or factor deficiencies in patient plasmas that test positive using ACTICLOT dPT and DVVtest.

Source plasmas for these products are processed in a manner consistent with established procedures to ensure the plasma is platelet-poor.¹

WARNING

The source material for these plasmas has been found to be non-reactive for Hepatitis B Surface Antigen (HBsAg), Hepatitis C Virus (HCV) and Human Immunodeficiency Virus Type 1 and Type 2 (HIV-1, HIV-2) using US Food and Drug Administration cleared methods. As no known test method can provide complete assurance that products derived from human blood will not transmit HBsAg, HCV, HIV-1, HIV-2 or other blood-borne pathogens, these plasmas should be handled as recommended for any potentially infectious human specimen.

REAGENT PREPARATION AND STABILITY

These lyophilised plasmas are stable until the date indicated on the vial label when stored at 2°C to 8°C. Reconstitute these plasmas with purified, deionised or distilled water as follows:

REF	816A	816N
Volume H ₂ O	0.5 mL	1.0 mL

After reconstitution, re-stopper the vial, swirl the contents gently and allow the vials to stand at room temperature for at least 15 minutes to ensure complete dissolution. Once reconstituted, these plasmas will remain stable for:

	8 hours
	2° - 8°C

Mix gently prior to use. Refer to the ACTICLOT dPT, *DVVtest* and *DVVconfirm* package inserts for complete assay procedures.

EXPECTED RESULTS/PERFORMANCE CHARACTERISTICS

LAtrol Abnormal Control and LAtrol Normal Control plasmas have been used in a quality control protocol for ACTICLOT dPT, *DVVtest* and *DVVconfirm* reagents. Assay values were generated using the ACL® 300R, BCT®, CA-7000, MLA® 900C, ST4 and STA Compact® coagulation analyzers. Instrument applications are available upon request. Clotting time results shown are for example only. The actual values will vary among laboratories and are dependent upon the method and instrumentation employed for clotting end-point detection. The protocol included multiple testing over several days. The overall coefficients of variation for the normal control and abnormal control are listed. The limits of quality control acceptability for these control plasmas should be established by each laboratory as part of a complete quality control protocol.

Table 1. LATrol™ Precision Data with ACTICLOT® dPT™

Coagulation Analyzer	Control	dPT Screening (sec)	Intra-Assay CV (%)	Inter-Assay CV (%)	dPT Confirmatory (sec)	Intra-Assay CV (%)	Inter-Assay CV (%)
ACL® 300R	816N	32.8	2.5	5.1	30.2	3.8	5.3
	816A	63.8	1.9	7.1	36.9	3.2	3.8
BCT®	816N	47.5	0.5	3.2	51.6	1.7	4.5
	816A	89.2	0.6	5.2	61.9	1.2	3.7
CA 7000	816N	44.3	0.9	ND	37.9	2.0	Nd
	816A	84.1	5.4	ND	48.6	3.3	Nd
MLA® 900C	816N	27.9	2.5	3.7	27.2	2.8	4.1
	816A	51.6	2.4	8.6	30.5	1.5	3.7
ST4	816N	38.1	0.7	ND	37.3	3.8	ND
	816A	60.6	6.5	ND	44.0	6.7	ND
STA Compact®	816N	40.2	0.8	3.4	39.7	0.9	4.3
	816A	77.9	1.1	7.2	46.0	1.0	4.8

ND – Not Determined

Table 2. LATrol™ Precision Data with DVVtest® and DVVconfirm®

Coagulation Analyzer	Control	DVVtest (sec)	Intra-Assay CV (%)	Inter-Assay CV (%)	DVVconfirm (sec)	Intra-Assay CV (%)	Inter-Assay CV (%)
ACL® 300R	816N	30.4	1.4	ND	31.3	0.6	ND
	816A	64.2	3.2	ND	35.1	1.4	ND
BCT®	816N	31.8	0.3	2.4	33.2	0.5	3.7
	816A	63.8	0.5	2.7	43.0	0.6	4.8
MLA® 900C	816N	31.0	0.9	2.9	31.7	1.3	2.3
	816A	63.6	1.4	3.9	34.5	0.9	3.2
ST4	816N	36.6	0.9	ND	34.4	1.0	ND
	816A	80.7	1.2	ND	41.1	1.2	ND
STA Compact®	816N	36.2	0.7	3.3	34.4	0.7	3.6
	816A	71.9	0.9	3.8	38.9	1.3	5.6

ND – Not Determined

REFERENCES

1. Brandt, J. T., Triplett, D. A., Alving, B., and Scharrer, I. Criteria for the Diagnosis of Lupus Anticoagulants, An Update. *Thrombosis and Haemostasis* 1995; 74(4): 1185-90.

ACL is a registered trademark of Instrumentation Laboratory, SpA

BCT is a registered trademark of Dade Behring Inc.

MLA is a registered trademark of Instrumentation Laboratory, SpA

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816A816N_A©ADI20110103
STA Compact is a registered trademark of Diagnostica Stago SA
ACL è un marchio commerciale di Instrumentation Laboratory, SpA
BCT è un marchio commerciale di Dade Behring Inc.
MLA è un marchio commerciale di Instrumentation Laboratory, SpA