# PLASMA FIBRINOGEN CONTROL LOW Ref SC070K

Abnormal Control Plasma for Fibrinogen coagulation assay quality control

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### ENGLISH

#### INTENDED USE:

Plasma Fibrinogen Control Low is a set of 6 vials of normal citrated human plasma proposed for the quality control of Fibrinogen clotting assays. The Fibrinogen concentration is measured using Fibrinogen assay (ref CK571K or CK572K or CK575K) from HYPHEN BioMed.

#### REAGENTS SUPPLIED:

6 vials of 1 ml of citrated human plasma, lyophilised.

**Note:** It contains an antibiotic as preservative. Each donor unit used for the preparation of this control plasma is a human plasma, which has been tested with registered methods for the presence of Hepatitis B Surface Antigen, Hepatitis C Antibodies (HVC) and antibodies to HIV 1 and 2 and was found negative. However, no test can completely exclude the presence of infectious agents. Any product of human origin, and more especially plasma, must be considered as being potentially infectious and must be handled with all the required cautions for this kind of material.

#### STORAGE CONDITIONS:

Unopened reagents, must be stored at 2–8 °C, in their original packaging box. They are then stable until the expiration date printed on the label.

**Note:** The stability studies at 30°C show that the reagent can be shipped at room temperature for a short period without damage.

#### Preparation and STABILITY OF REAGENTS:

Reconstitute each vial with exactly 1 mL of distilled water. Shake thoroughly until complete dissolution of the content (vortex). Incubate at room temperature (18–25°C) for 30 min, while shaking the vial from time to time. Homogenise the content before each use. Stability of the control plasma, kept in its original vial:

- 8 hours at room temperature (18–25°C)
- or 24 hours at 2–8°C.
- Do not freeze.

#### Cautions:

- In order to improve stability, reagents must be closed with their original screw caps following each use.
- Reagents must be handled with care, in order to avoid any contamination during use.
- It is recommended to homogenize each vial before use, in order to have a good reproducibility, all the time.

### APPLICATIONS:

Plasma Fibrinogen Control Low can be used as an abnormal control plasma for Fibrinogen clotting assays. The Fibrinogen concentration reported is accurately determined by reference to an internal reference standard, initially validated by reference to the corresponding NIBSC (National Institute for Biological Standards and Controls, UK) International Standard. The following table shows the usual ranges expected for the Plasma Fibrinogen Control Low.

Parameter	Method	Usual range for the target value
Fibrinogen	Clotting	1.2 – 1.7 g/L

### PERFORMANCE CHARACTERISTICS:

The following values, obtained for one lot, are provided as an example only.

Parameter	[C]	Acceptance range
Fibrinogen	1.47	1.17 - 1.77

The concentration and acceptance range may present variations from lot to lot, but it is exactly measured for each lot and reported on the **flyer** provided within the kit.

When the Control is used as quality control plasma for the assay of Fibrinogen by clotting method, the values obtained must be within the acceptance ranges reported for the lot used, in order to validate the test series. Should the value be out of these ranges, the results for the corresponding series must be considered as invalid. It is then recommended to rerun the series and to check all the assay parameters.

**Note:** The results can slightly vary according to the reagent and the instrument used. The target value and acceptance range must consequently be confirmed and adjusted, if necessary, for each new lot of control, in the laboratory working conditions.

### CAUTIONS:

- Like all lyophilised plasmas, the plasmas from the Plasma Fibrinogen Control Low are more or less cloudy after reconstitution. This is due essentially to the lipids that, after lyophilisation, become less soluble and can form a light deposit.
- If necessary, let each vial 10 minutes at room temperature and shake gently before use in order to homogenise the content.
- Reagents must be handled with care, in order to avoid any contamination or activation during use. Any plasma containing a coagulum or contamination must be rejected.

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