Calcium Chloride 0.025M (M/40)

REF AR001K, 12x25 mL

REF AR001B, 3x15 mL

REF AR001L, 12x15 mL

Calcium Chloride 0.025M (M/40)

Chlorure de Calcium 0,025M (M/40)

English, last revision: 01-2021

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INTENDED USE:

Calcium Chloride solution for coagulation assays, such as the Activated Partial Thromboplastin Time (aPTT) or for the assays of factors of the intrinsic coagulation pathway.
Ready to use.

REAGENTS:

Calcium Chloride Solution. Contains Calcium Chloride (0.025M) and small amounts of Sodium azide (0.9 g/L) as preservative.

REF AR001K → 12 vials of 25 mL

REF AR001B → 3 vials of 15 mL

REF AR001L → 12 vials of 15 mL

WARNING AND PRECAUTIONS:

- In contact with lead or copper pipes, sodium azide can generate explosive compounds.
- Waste should be disposed of in accordance with applicable local regulations.
- Use only the reagents from the same batch of kits.
- Aging studies show that the reagents can be shipped at room temperature without degradation.
- This device of *in vitro* diagnostic use is intended for professional use in the laboratory.

REAGENT PREPARATION:

Reagent is ready to use, homogenize and load it on the analyzer following application guide instruction.

For another use, allow to stabilize for 30 minutes at room temperature (18-25°C), homogenize before use.

STORAGE AND STABILITY:

Unopened reagents should be stored at 2-8°C in their original packaging. Under these conditions, they can be used until the expiry date printed on the kit.

Reagent stability after opening, free from any contamination or evaporation, and stored closed, is of:

- 8 weeks at 2-8°C.
- 8 weeks at room temperature (18-25°C).
- Stability on board of the analyzer: see the specific application.

LIMITATIONS:

- Any reagent presenting an unusual appearance or showing signs of contamination must be rejected.
- Refer to the specific associated kit inserts for instructions concerning the test to be performed.

SYMBOLS:

Symbols used and signs listed in the ISO 15223-1 standard, see Symbol definitions document.

Changes compared to the previous version.

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